











HOW TO SET UP TEST PATIENTS IN ARIA FOR REGIMEN VALIDATION

As part of the regimen validation process, two test patients should be set-up for each regimen to be tested, one for the pharmacist validation and one for clinician validation which is then used for the nurse validation.

Adult Regimens

<p>1. New Patient</p> <p>Once logged in, Click New under the Patient tab</p>																				
<p>2. PHARMACIST TEST PATIENT</p> <p>New → General Tab</p> <ul style="list-style-type: none"> ● Patient name: use K&M SACT protocol number, e.g. BRA-002, as follows: <ul style="list-style-type: none"> ○ Last name (upper case): XXBRA ○ First name: 002 ● D.O.B.: enter any date but ensure patient will be over 18 years ● Sex: Enter either Male or Female ● Error/test patient: radio button must be checked 	<p>CONSULTANT TEST PATIENT</p> <p>New → General Tab</p> <ul style="list-style-type: none"> ● Patient name: use K&M SACT protocol number, e.g. BRA-002, as follows: <ul style="list-style-type: none"> ○ Last name (upper case): XXTEST002 ○ First name: BRA ● D.O.B.: enter any date but ensure patient will be over 18 years ● Sex: Enter either Male or Female ● Error/test patient: radio button must be checked 																			
<p>3. New → IDs Tab</p> <p>R number: BRA-002</p> <p>If there are multiple test patients for the same regimen, version number can be added.</p> <p>Click OK</p>	<p>New → IDs Tab</p> <p>R number: BRA002</p> <p>If there are multiple test patients for the same regimen, version number can be added.</p> <p>Click OK</p>																			
<p>4.</p> <p>Open the patient you've created by double clicking on their entry.</p>	 <table border="1"> <thead> <tr> <th>Patient</th> <th>R Number</th> <th>Date of Birth</th> </tr> </thead> <tbody> <tr> <td>xxTest002, Gyn</td> <td>gyn002</td> <td>May 26, 1963</td> </tr> <tr> <td>XXTEST002, BRA</td> <td>BRA002v6</td> <td>Aug 27, 1965</td> </tr> <tr style="background-color: #e0e0e0;"> <td>XXTEST002, BRAv2</td> <td>BRA002v2</td> <td>Apr 20, 1956</td> </tr> <tr> <td>XXTEST002, HAEMCML</td> <td>HAEMCML002</td> <td>Jan 1, 1950</td> </tr> <tr> <td>XXTest002+023, GYN</td> <td>GYN002+023</td> <td>Jan 1, 1965</td> </tr> </tbody> </table>		Patient	R Number	Date of Birth	xxTest002, Gyn	gyn002	May 26, 1963	XXTEST002, BRA	BRA002v6	Aug 27, 1965	XXTEST002, BRAv2	BRA002v2	Apr 20, 1956	XXTEST002, HAEMCML	HAEMCML002	Jan 1, 1950	XXTest002+023, GYN	GYN002+023	Jan 1, 1965
Patient	R Number	Date of Birth																		
xxTest002, Gyn	gyn002	May 26, 1963																		
XXTEST002, BRA	BRA002v6	Aug 27, 1965																		
XXTEST002, BRAv2	BRA002v2	Apr 20, 1956																		
XXTEST002, HAEMCML	HAEMCML002	Jan 1, 1950																		
XXTest002+023, GYN	GYN002+023	Jan 1, 1965																		
<p>5. History Window</p> <ul style="list-style-type: none"> ● Allergies Tab: Check 'No Known Allergies' ● Medications Tab: Check 'No External Medications' 																				
<p>6. Vital Signs Window</p> <ul style="list-style-type: none"> ● Large patient <ul style="list-style-type: none"> ○ Height 190cm ○ Weight 115kg <p>Click All Approve</p>																				

Paediatric Regimens

<p>1. New Patient</p> <p>Once logged in, Click New under Patient</p>																										
<p>2. PHARMACIST TEST PATIENT</p> <p>New → General Tab</p> <ul style="list-style-type: none"> • Patient name: use K&M SACT protocol number, e.g. BRA-002, as follows: <ul style="list-style-type: none"> ○ Last name (upper case): XXBRA ○ First name: 002 • D.O.B.: enter any date but ensure patient will be under 18 years • Sex: Enter either Male or Female • Error/test patient: radio button must be checked 	<p>CONSULTANT TEST PATIENT</p> <p>New → General Tab</p> <ul style="list-style-type: none"> • Patient name: use K&M SACT protocol number, e.g. BRA-002, as follows: <ul style="list-style-type: none"> ○ Last name (upper case): XXTEST002 ○ First name: BRA • D.O.B.: enter any date but ensure patient will be under 18 years • Sex: Enter either Male or Female • Error/test patient: radio button must be checked 																									
<p>3. New → IDs Tab</p> <p>R number: BRA-002</p> <p>If there are multiple test patients for the same regimen, version number can be added.</p> <p>Click OK</p>	<p>New → IDs Tab</p> <p>R number: BRA002</p> <p>If there are multiple test patients for the same regimen, version number can be added.</p> <p>Click OK</p>																									
<p>4.</p> <p>Open the patient you've created by double clicking on their entry.</p>	<table border="1"> <thead> <tr> <th>Patient</th> <th>R Number</th> <th>Date of Birth</th> <th></th> </tr> </thead> <tbody> <tr> <td>xxTest002, Gyn</td> <td>gyn002</td> <td>May 26, 1963</td> <td></td> </tr> <tr> <td>XXTEST002, BRA</td> <td>BRA002v6</td> <td>Aug 27, 1965</td> <td></td> </tr> <tr style="border: 2px solid orange;"> <td>XXTEST002, BRAv2</td> <td>BRA002v2</td> <td>Apr 20, 1956</td> <td></td> </tr> <tr> <td>XXTEST002, HAEMCML</td> <td>HAEMCML002</td> <td>Jan 1, 1950</td> <td></td> </tr> <tr> <td>XXTest002+023, GYN</td> <td>GYN002+023</td> <td>Jan 1, 1965</td> <td></td> </tr> </tbody> </table>		Patient	R Number	Date of Birth		xxTest002, Gyn	gyn002	May 26, 1963		XXTEST002, BRA	BRA002v6	Aug 27, 1965		XXTEST002, BRAv2	BRA002v2	Apr 20, 1956		XXTEST002, HAEMCML	HAEMCML002	Jan 1, 1950		XXTest002+023, GYN	GYN002+023	Jan 1, 1965	
Patient	R Number	Date of Birth																								
xxTest002, Gyn	gyn002	May 26, 1963																								
XXTEST002, BRA	BRA002v6	Aug 27, 1965																								
XXTEST002, BRAv2	BRA002v2	Apr 20, 1956																								
XXTEST002, HAEMCML	HAEMCML002	Jan 1, 1950																								
XXTest002+023, GYN	GYN002+023	Jan 1, 1965																								
<p>5. History Window</p> <ul style="list-style-type: none"> • Allergies Tab: Check 'No Known Allergies' • Medications Tab: Check 'No External Medications' 																										
<p>6. Vital Signs Window</p> <ul style="list-style-type: none"> • Large patient <ul style="list-style-type: none"> ○ Height 150cm ○ Weight 50kg <p>Click All Approve</p>	